BEAUTIFLY

Therm Body Revive

Golden

series

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[patrz str. 3, f7]

| | |
|-------------------|-----------------|
| Input voltage: | 5V/1A |
| Size: | 130*86*58 |
| Libra: | 316g |
| Charging time: | 4 hours |
| Working time: | 15 minutes |
| Battery type: | lithium battery |
| Battery capacity: | 2500mAh |

Included [see page 2. f1]

1. Device
2. User manual
3. Warranty card
4. USB charging cable
5. Vacuum massage pads

Device description [see page 2. f2]

1. Charging port
2. Particulate filter
3. Compress

Panel display [see page 2. f3]

1. Numerical indicator
2. Battery indicator
3. Heating indicator
4. Pressure indicator
5. Heart rate indicator

Operation panel [see page 2. f4]

1. On/Off
2. Pressure button
3. Heating button
4. Intensity + button
5. Intensity button -

Control panel description [see page 2. f5]

1. On/Off - press and hold for 3 seconds to turn on the device. A short press changes the vacuum intensity, from 1 to 5 levels.
2. Pressure Hold/Pressure Release - the first short press holds pressure, the second short press releases pressure, and the third short press creates a vacuum.
3. Heating - short press of the heating button changes the temperature from level 1 to level 3. (40°C, 45°C, 50°C)
4. Negative pressure button - press the button to turn on, press again to turn off.
5. Pressure intensity button - press the button to turn on, short press changes the intensity, between 1 and 5 levels of intensity.

User manual

1. Press the "on/off" button for 2 seconds to turn on the device, then short press the button to select the intensity level that suits you, from level 1 to level 5.
2. If you need to turn on the heating function, press the "Heating" button and select a level from 1 to 3 levels that is suitable for you.
3. After use, short press "Pressure/Release" and release the pressure to turn off the device.

Instructions - massage with accessories [see page 3. f6]

1. Press "on/off" for 2 seconds to turn on the device and select the appropriate suction cup.
2. Fit the suction cup to the device, and then select the appropriate massage intensity.
3. After use, pull back the pressure relief valve on the device and remove it.
4. Press "on/off" for 2 seconds to turn off the device.

Storage

- Keep the device dry and clean, out of the reach of children.
- Store away from high temperatures.
- Store the device in a cool, dry place.

Cleaning

- Keep the device clean. For daily maintenance, clean it with a neutral detergent and a soft cloth. Do not clean the device with organic solvents or other chemicals. Do not rinse it with water.
- Avoid high pressure washing and scratching the surface of the product with hard objects.

Precautions and safety

- The operating environment of the device should be well ventilated, the temperature should be kept at 5°C~40°C, and the humidity should be 30%~85%.
- Keep the device away from electrical appliances, sources of high-frequency electrical radiation, etc., to avoid interference during use.
- Under no circumstances should you disassemble, repair or modify the device.
- When using for the first time, time should be controlled (5-15 min). Moreover, pay attention to the frequency. Excessive use may cause discomfort. Stop using the device if you feel uncomfortable.
- In case of incorrect operation of the device, stop using it immediately and contact the service.

Warnings and contraindications

- If the device has not been used for a long time, wipe it down and turn it on to check if it works properly. Then charge it. If any irregularities occur, contact the service.
- Persons with implanted electronic devices, such as pacemakers, must not use the device.
- Using the device near a short wave or microwave processing device (eg 1M) may cause the output of the device to become unstable.
- Do not use this device during pregnancy and menstruation.
- Keep the device and its accessories out of the reach of children.
- Do not immerse the device in any liquid and do not wipe the product with a wet cloth.
- Please turn off the power before charging the device.
- Connect to USB for charging. The device cannot be used while it is charging.
- Improper use or adjustment of the device may void the service warranty.
- Persons with wounds, fractures, scars, cancer, varicose veins and local poor skin elasticity, must not use the device.
- Persons with serious heart, liver, kidney diseases or high fever, must not use the device.

More information about the product is available on the manufacturer's website at: [see page 3, f7]

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EN: Electrical and electronic devices, as well as batteries and accumulators, must not be mixed with other household waste. This is indicated by the crossed-out trash can symbol on the packaging. Used electrical and electronic devices, batteries, and accumulators must be taken to local selective waste collection points or to the retailer, and more information can be found in your municipality. Electrical and electronic equipment, batteries, and accumulators may contain hazardous substances, mixtures, and components harmful to the environment and human health. Improper handling or damage to them can lead to health hazards or environmental pollution during further disposal or recycling.

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EN: The manufacturer hereby declares, at its sole responsibility, that the Therm Body Revive Golden product complies with the Directive 2014/30/UE of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as with the Directive 2011/65/UE of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/UE with regard to the list of restricted substances. The full text of the declaration of conformity is available at the following website:

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